




F-System

Bonegrafting Screws




Table of content

Introduction	3
Features	4
Product information	5


F-series: Bonegrafting System

INTRODUCTION

Bone grafting can repair implant sites with inadequate bone structure due to previous extractions, gum disease or injuries. The bone is either obtained from the jaw, hip or cranial. Sinus bone grafts are also performed to replace bone in the posterior upper jaw.



1. Simple and efficient
2. Compatibility
3. Small in all dimensions



INDICATIONS

PRINCIPLES



F-series: Bonegrafting System

FEATURES



SCREWHEAD



- Easy pickup with square screwdriver IT0128 (or IT0127)
- Easy removal with hexagon screwdriver IT0129 (or IT0123)



SCREW THREAD




- An excellent grip on the bone
- The thread prevents stripping through the bone
- Excellent retention

TIP



- No damage to the soft tissue

BOX



The screws are nicely arranged in a sterilizable box.

F-series: Bonegrafting System

PRODUCT INFORMATION



IT0123




IT0128




IT0127



IT0129



IT3002




IT3007

Cat. No.	Description
IT0123	Square Screwdriver 1.2mm - LONG
IT0128	Screwdriver Internal Hex. long, titanium
IT0127	Square Screwdriver 1.2mm - SHORT
IT0129	Titanium screwdriver 1.2 mm
IT3002	Organiser F-series 4 pc
IT3007	Organiser F-series 3 pc



Cat. No.	Description
F1203	Titanium miniscrew 1.2 x 3mm
F1204	Titanium miniscrew 1.2 x 4 mm
F1205	Titanium miniscrew 1.2 x 5 mm
F1206	Titanium miniscrew 1.2 x 6 mm
F1208	Titanium miniscrew 1.2 x 8 mm
F1210	Titanium miniscrew 1.2 x 10 mm
F1212	Titanium miniscrew 1.2 x 12 mm
F1214	Titanium miniscrew 1.2 x 14 mm



INSTRUMENTS

SCREWS



F-series: Bonegrafting System

PRODUCT INFORMATION



TWIST DRILLS




Cat. No.	Description
IT105007	Twist drill, Ø 1.0 mm, L=50 mm, stop 7 mm
IT105005	Twist drill, Ø 1.0 mm, L=50 mm, stop 5 mm



Cat. No.	Description
IT105014	Twist drill Ø 1.0mm, L=50 mm, stop 14 mm

BOX CONTENT




IT3007


IT3007	Cat. No.	Description	Max. quantity
➤	IT0127	Square Screwdriver 1.2mm - SHORT	2
➤	IT0129	Titanium screwdriver 1.2 mm	1
➤	IT105014	Twist drill Ø 1.0mm, L=50 mm, stop 14 mm	2
➤	F1203	Titanium miniscrew 1.2 x 3mm	10
➤	F1204	Titanium miniscrew 1.2 x 4 mm	10
➤	F1205	Titanium miniscrew 1.2 x 5 mm	10
➤	F1206	Titanium miniscrew 1.2 x 6 mm	10
➤	F1208	Titanium miniscrew 1.2 x 8 mm	10
➤	F1210	Titanium miniscrew 1.2 x 10 mm	10
➤	F1212	Titanium miniscrew 1.2 x 12 mm	10
➤	F1214	Titanium miniscrew 1.2 x 14 mm	10

F-series: Bonegrafting System

PRODUCT INFORMATION




IT3002



BOX CONTENT

IT3002	Cat. No.	Description	Max. quantity
➤	IT0123	Square Screwdriver 1.2mm - LONG	1
➤	IT0127	Square Screwdriver 1.2mm - SHORT	1
➤	IT0128	Screwdriver Internal Hex. long, titanium	2
➤	IT105014	Twist drill Ø 1.0mm, L=50 mm, stop 14 mm	2
➤	F1203	Titanium miniscrew 1.2 x 3mm	10
➤	F1204	Titanium miniscrew 1.2 x 4 mm	10
➤	F1205	Titanium miniscrew 1.2 x 5 mm	10
➤	F1206	Titanium miniscrew 1.2 x 6 mm	10
➤	F1208	Titanium miniscrew 1.2 x 8 mm	10
➤	F1210	Titanium miniscrew 1.2 x 10 mm	10
➤	F1212	Titanium miniscrew 1.2 x 12 mm	10
➤	F1214	Titanium miniscrew 1.2 x 14 mm	10





Headquarters:

Titamed
Mechelsesteenweg 315 G
2550 Kontich
Belgium

tel: +32 (0)3 289 51 56
email: info@titamed.be
website: www.titamed.com

SP-605 Rev A

