

Modules can be changed or removed without removing the plates and the screws.

Very useful distractor for pathological small upper jaws and clefts (module 8, 10 en 12mm).


INOX CONTAINER



SMILE 3 DISTRACTOR MODULE



ABUTMENT PLATES




SELF-DRILLING SCREWS



BLOCKING SCREW



INSTRUMENTS



Cat. No.	Description
IT.3003	Inox Container for Smile Distractor
IT.1011	Drill with round tip

Cat. No.	Length of distractor (mm)	Total length of expansion
DS.008(A)	8	8+6
DS.010(A)	10	10+8
DS.012(A)	12	12+10
DS.014(A)	14	14+12
DS.016(A)	16	16+14
DS.018(A)	18	18+16
DS.020(A)	20	20+18

Cat. No.	Description
D001	Abutmentplates

Cat. No.	Description
DSSA06	6 mm, Ø 2.3 mm
DSSA07	7 mm, Ø 2.3 mm

Cat. No.	Description
D.002	Ø 2.0 mm

Cat. No.	Description
IT.0315	Osteosynthesis Screwdriver
IT.0311	Form Stamp Right
IT.0312	Form Stamp Left
IT.0314	Distractor Key Instrument
IT.0313	Distractor Holder
IT.0128	Distractor Blocking Screwdriver
IT.0333	Patient Key