

Modules can be changed or removed without removing the plates and the screws.

Very useful distractor for pathological small upper jaws and clefts (module 8, 10 en 12mm).


### INOX CONTAINER



### SMILE 3 DISTRACTOR MODULE



### ABUTMENT PLATES




### SELF-DRILLING SCREWS



### BLOCKING SCREW



### INSTRUMENTS



Cat. No.	Description
IT.3003	Inox Container for Smile Distractor
IT.1011	Drill with round tip

Cat. No.	Length of distractor (mm)	Total length of expansion
DS.008(A)	8	8+6
DS.010(A)	10	10+8
DS.012(A)	12	12+10
DS.014(A)	14	14+12
DS.016(A)	16	16+14
DS.018(A)	18	18+16
DS.020(A)	20	20+18

Cat. No.	Description
D001	Abutmentplates

Cat. No.	Description
DSSA06	6 mm, Ø 2.3 mm
DSSA07	7 mm, Ø 2.3 mm

Cat. No.	Description
D.002	Ø 2.0 mm

Cat. No.	Description
IT.0315	Osteosynthesis Screwdriver
IT.0311	Form Stamp Right
IT.0312	Form Stamp Left
IT.0314	Distractor Key Instrument
IT.0313	Distractor Holder
IT.0128	Distractor Blocking Screwdriver
IT.0333	Patient Key